




Product Datasheet

CD133 antibody orb99113

Description	CD133 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD133. Please contact us for the exact immunogen sequence . The peptide is available as orb374954.
Target	CD133. Human;Mouse;Rat. No cross reactivity with other proteins.
Form/Appearance	10 mM PBS, 0.02% sodium azide, 10 mg/ml rAlbumin
Concentration	- 100 µg (in 100 µl): 1 mg/ml - 200 µg (in 200 µl): 1 mg/ml
Storage	After reconstitution store at 4°C. for up to two weeks. For long term storage aliquot and store at -20°C. Avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	95 kDa
Uniprot ID	O54990 , O43490
NCBI	1 , 001145847 , 001139319 , 11
Entrez	8842
Dilution Range	IHC-P:1:200, WB: 1:500, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.



IF analysis of mouse skin tissue using CD133 antibody (Dilution of primary antibody 1:200)



IHC-P image of human endometrial cancer tissue using anti-CD206 (dilution of primary antibody at 1:100)