




Product Datasheet

TIMP1 antibody orb26245

Description	TIMP1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Metalloproteinase inhibitor 1 protein (24-207AA)
Target	TIMP1
Form/Appearance	Liquid
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01033
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human breast cancer using TIMP1 antibody



Western blot analysis of Mouse heart using TIMP1 antibody