



## Product Datasheet

### ApoE Antibody orb69429

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to ApoE                                 |
| <b>Species/Host</b>        | Mouse                                                             |
| <b>Reactivity</b>          | Human                                                             |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Tested Applications</b> | ELISA, FC, IHC                                                    |
| <b>Immunogen</b>           | Purified recombinant fragment of human ApoE expressed in E. Coli. |
| <b>Concentration</b>       | 1mg/ml                                                            |
| <b>Preservatives</b>       | Ascitic fluid containing 0.03% sodium azide.                      |
| <b>Storage</b>             | 4°C; -20°C for long term storage. Avoid freeze/thaw cycles.       |
| <b>Note</b>                | For research use only.                                            |
| <b>Isotype</b>             | Mouse IgG1                                                        |
| <b>Clonality</b>           | Monoclonal                                                        |
| <b>MW</b>                  | 36kDa                                                             |
| <b>Clone ID</b>            | 1H4A1                                                             |
| <b>Entrez</b>              | <a href="#">348</a>                                               |
| <b>Dilution Range</b>      | IHC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |



Western blot analysis of HEK293 (Lane1) and ApoE (AA: 20-267)-hlgGfc transfected HEK293 (Lane2) cell lysate using ApoE antibody

