

Product Datasheet

BrdU antibody (FITC) orb461471

Description	BrdU antibody (FITC)
Species/Host	Mouse
Reactivity	Human, Rat
Predicted Reactivity	Other
Conjugation	FITC
Tested Applications	IF
Immunogen	KLH conjugated Brdu
Target	BrdU(Proliferation Marker)
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Monoclonal
MW	0.3071kDa
Clone ID	1A7
Dilution Range	IF=1:100-500
Expiration Date	12 months from date of receipt.