



Product Datasheet

TIMP1 antibody orb195994


Description	Rabbit polyclonal antibody to TIMP1
Species/Host	Rabbit
Reactivity	Bovine, Canine, Guinea pig, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human TIMP1. Please contact us for the exact immunogen sequence. The peptide is available as orb374835 .
Target	TIMP1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml - 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	P01033 , P12032 , P30120
NCBI	2 , 003254 , 21 , 1 , 003245 , 11
Entrez	7076
Dilution Range	WB: 1:100-800, IHC-P: 1:400, IF/ICC: 1:400
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat stomach tissue using TIMP1 antibody (2.5ug/ml)



Immunohistochemical staining of paraffin embedded mouse skin tissue using TIMP1 antibody (primary antibody at 1:100)



Immunohistochemical staining of rat prostate tissue using TIMP1 antibody (2.5ug/ml)