



Product Datasheet

APP antibody orb216018

Description	Rabbit polyclonal antibody to Beta Amyloid (Phospho-T743)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, IF, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T743 of human Beta Amyloid protein. The exact sequence is proprietary.
Target	APP
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P12023 , P08592 , P05067
Entrez	11820 , 351 , 54226
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF7 colchicine-treated (Lane 1), mouse liver (Lane 2), PC12 (Lane 3) whole cell lysates using APP antibody



Immunohistochemical analysis of formalin-fixed and paraffin-embedded human breast cancer tissue using APP antibody (dilution at 1:200)