

Product Datasheet

Somatostatin Receptor 2 antibody (FITC) orb222574

Description	Rabbit polyclonal antibody against Somatostatin Receptor 2 conjugated to FITC.
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Bovine, Canine, Equine, Guinea pig, Mouse, Porcine, Rabbit, Rat
Conjugation	FITC
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human Somatostatin Receptor 2
Target	SSTR2
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	41kDa
Uniprot ID	P30874
Dilution Range	IF=1:100-500
Expiration Date	12 months from date of receipt.