

## Product Datasheet

### Recombinant Human FGF10 Protein orb527080

|                  |                                                                                                                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity       | Human                                                                                                                                                                                                                  |
| Form/Appearance  | Lyophilized                                                                                                                                                                                                            |
| Storage          | Lyophilized recombinant human KGF-2/FGF10 remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, KGF-2/FGF10 remains stable up to 2 weeks at 4°C or up to 3 months at -20°C.                |
| Note             | For research use only.                                                                                                                                                                                                 |
| Protein Sequence | QALGQ DMVSP EATNS SSSSF SSPSS AGRHV RSYNH LQGDV<br>RWRKL FSFTK YFLKI EKNGK VSGTK KENCY YSILE ITSVE IGVVA<br>VKAIN SNYYL AMNKK GKLYG SKEFN NDCKL KERIE ENGYN TYASF<br>NWQHN GRQMY VALNG KGAPR RGQKT RRKNT SAHFL PMVVH S |
| MW               | 19.4KD                                                                                                                                                                                                                 |
| Source           | E. coli Gln38-Ser208                                                                                                                                                                                                   |
| Uniprot ID       | <a href="#">O15520</a>                                                                                                                                                                                                 |
| Expiration Date  | 6 months from date of receipt.                                                                                                                                                                                         |