

Product Datasheet

Dehydroepiandrosterone-BSA conjugate orb527316

Description	Dehydroepiandrosterone-BSA conjugate
Tested Applications	ICA, ELISA, GICA
Hazard Information	For research use only
Form/Appearance	Dry powder
Preservatives	0.01mol/L, pH 7.4 PBS
Storage	Short-term store at 4?. Long-term store at -20?. Avoid repeated freeze.
Note	For research use only.
CAS Number	53-43-0
Source	Bovine serum albumin (BSA)
Expiration Date	12 months from date of receipt.
Application Notes	Other Small Molecule Antigen