

Product Datasheet

Retinoic acid- HSA conjugate orb527386

Description	Retinoic acid- HSA conjugate
Tested Applications	ICA, ELISA
Hazard Information	For research use only
Form/Appearance	Dry powder
Preservatives	0.01mol/L, pH 7.4 PBS
Storage	Short-term store at 4?. Long-term store at -20?. Avoid repeated freeze.
Note	For research use only.
CAS Number	302-79-4
Source	Human serum albumin (HSA)
Expiration Date	12 months from date of receipt.
Application Notes	Other Small Molecule Antigen