



Product Datasheet

TIMP1 antibody orb542013


Description	Rabbit polyclonal antibody to TIMP1
Reactivity	Human, Mouse, Rat
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human TIMP1(NP_003245.1).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	26KDa
Uniprot ID	P01033
NCBI	7076
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cells lysates using TIMP1 antibody



Western blot analysis of extracts of mouse lung lysates using TIMP1 antibody



Western blot analysis of extracts of various cell lines lysates using TIMP1 antibody