

Product Datasheet

Opti-Protein Ultra Marker orb545606

Description	Opti-Protein Ultra Marker
Hazard Information	Non-Toxic
Target	Opti-Protein Ultra Marker
Form/Appearance	Lyophilized Powder
Storage	Store 4°C for up to 3 months or -20°C for up to 24 months.
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	• Monitoring protein migration during SDS-PAGE • Monitoring efficiency of protein transfer onto PVDF or nitrocellulose for Western Blotting • Sizing of proteins on SDS-PAGE and Western Blots