


Product Datasheet

Mouse TNF Alpha/Tumor necrosis factor ELISA Kit orb50113

Description	Mouse TNF Alpha/Tumor necrosis factor ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Sensitivity	
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: L80-L235
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only.
MW	25896 MW
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	P06804
Expiration Date	Please enquire.



Mouse TNF alpha ELISA
Kit