

---

## Product Datasheet

### GR Kit (Total RNA Isolation Reagent) orb180503

#### Description

The GR Kit provides an efficient 3-step method to isolate the total RNA from the tissue, cultured animal and bacterial cells, blood, and serum. This unique reagent system ensures the total RNA with a high yield and good quality from samples of unlimited size. If a larger sample is required, the reagent volume can be scaled proportionately, making this reagent not only very user-friendly but also highly versatile. The RNA phenol extraction is not required, and the entire procedure can be completed in 60 minutes. The total RNA is ready for use in RT-PCR, Northern Blotting, cDNA Synthesis and Mapping.

#### Note

For research use only.

#### Expiration Date

12 months from date of receipt.

#### Application Notes

Nucleic Acid Purification System