

Product Datasheet

Mouse CD120b ELISA Kit orb219607

Description	Sandwich ELISA for quantitative detection of Mouse CD120b concentrations.
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Sensitivity	1 pg/ml
Tested Applications	ELISA
Target	TNFRSF1B
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only.
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P25119
Entrez	21938
Expiration Date	Please enquire.
Application Notes	Cell Culture Supernatants, SerumExperiment Notes: Colorimetric