

## Product Datasheet

### Rat VEGFA ELISA Kit orb391179

|                     |                                                                                                          |
|---------------------|----------------------------------------------------------------------------------------------------------|
| Description         | Sandwich ELISA for quantitative detection of rat VEGFA concentrations.                                   |
| Reactivity          | Rat                                                                                                      |
| Range               | 15.6 pg/ml - 1000 pg/ml                                                                                  |
| Sensitivity         | 4 pg/ml                                                                                                  |
| Tested Applications | ELISA                                                                                                    |
| Target              | VEGFA                                                                                                    |
| Storage             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.            |
| Note                | For research use only.                                                                                   |
| Sample Types        | Cell Culture Supernates, Serum, Plasma, Tissue Homogenates                                               |
| Uniprot ID          | <a href="#">P16612</a>                                                                                   |
| Entrez              | <a href="#">83785</a>                                                                                    |
| Expiration Date     | Please enquire.                                                                                          |
| Application Notes   | Detection Method: ColorimetricSample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates |