


Product Datasheet

Complement C3d Fragment antibody orb608034

Description	Complement C3d Fragment antibody
Species/Host	Rabbit
Reactivity	Mouse
Predicted Reactivity	Human, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Complement C3d fragment (991-1090/1663 aa)
Target	Complement C3d fragment
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	33/109/181 kDa
Uniprot ID	P01024
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



WB analysis of mouse spleen lysate at 40 ug using Complement C3d Fragment antibody