



Product Datasheet

RNF105 antibody orb608057

Description	RNF105 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Predicted Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human RNF105 (251-350/2025 aa)
Target	RNF105
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	230 kDa
Uniprot ID	P53804
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of rat brain using RNF105 Polyclonal Antibody at 1:200 dilution.



IHC-P analysis of mouse brain using RNF105 Polyclonal Antibody at 1:200 dilution.