

Product Datasheet

IFI78 antibody orb608097

Description	IFI78 antibody
Species/Host	Rabbit
Predicted Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from mouse IFI78/MX1 (261-360/631 aa)
Target	MX1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	72 kDa
Uniprot ID	P20591
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.