

Product Datasheet

TNF-alpha antibody (HRP) orb608135

Description	TNF-alpha antibody (HRP)
Species/Host	Mouse
Predicted Reactivity	Insect, Monkey, Mouse, Rat, Yeast
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Recombinant human TNF alpha 77-233aa
Target	TNF alpha
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG1
Clonality	Monoclonal
Clone ID	3A6
Uniprot ID	P01375
Dilution Range	WB=1:1000-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.