


Product Datasheet

Bcl-xL (phospho-Ser62) antibody orb608144

Description	Bcl-xL (phospho-Ser62) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Predicted Reactivity	Bovine, Canine, Equine, Mouse, Porcine, Sheep
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Bcl-xL around the phosphorylation site of Ser62 (55-200/233 aa)
Target	phospho-Bcl-xL (Ser62)
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	26 kDa
Uniprot ID	Q07817
Dilution Range	IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



ICC analysis of Hela cell using phospho-Bcl-xL (Ser62) polyclonal Antibody at 1:100 dilution; DAPI (blue) was used to stain the cell nuclei.