

Product Datasheet

NCF4 antibody orb608148

Description	NCF4 antibody
Species/Host	Rabbit
Reactivity	Mouse
Predicted Reactivity	Bovine, Equine, Human, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human NCF4 (1-100/339 aa)
Target	NCF4
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q15080
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.