

Product Datasheet

PGK1 antibody orb608149

Description	PGK1 antibody
Species/Host	Rabbit
Reactivity	Mouse
Predicted Reactivity	Bovine, Equine, Human, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human PGK1 (241-340/417 aa)
Target	PGK1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	P00558
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.