

Product Datasheet

Methyltransferase Dilution Buffer III [Out of stock] orb608789

Description	Methyltransferase Dilution Buffer III [Out of stock]
Tested Applications	Kinase Assay
Hazard Information	Non-Toxic
Target	Methyltransferase Dilution Buffer III
Preservatives	Frozen solution of buffer containing: 50mM Tris-HCl, pH 7.5, 1mM EDTA, 5% Glycerol and 10ng/ul rAlbumin.
Note	For research use only.
Source	Solution of (1X) buffer used for diluting enzyme DNA methyltransferases.
Expiration Date	12 months from date of receipt.