

Product Datasheet

GSTO1 antibody orb610930

Description	Rabbit polyclonal antibody to GSTO1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human GSTO1
Target	GSTO1
Form/Appearance	Liquid
Concentration	Batch dependent within range: 100 ul at 0.5 - 1 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	P78417
Expiration Date	12 months from date of receipt.