

Product Datasheet

mCherry antibody orb611194

Description	Chicken polyclonal antibody to mCherry (Cherry fluorescent protein).
Species/Host	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Affinity purified recombinant Red Fluorescent Protein (mCherry).
Target	mCherry
Form/Appearance	Polyclonal antibody supplied as a 100 µl (2 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. Affinity purified using the immunogen immobilized on a solid support; Sample purity > 95% (SDS-PAGE).
Concentration	2 mg/ml
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only.
Isotype	IgY
Clonality	Polyclonal
Dilution Range	WB:1:10,000-1:30,000, IF:1:200-1:1,000, ELISA:1:20,000-1:50,000
Expiration Date	12 months from date of receipt.
Application Notes	The antibody solution should be gently mixed before use.