

Product Datasheet

Human QPCT ELISA Kit orb610884

Description Human QPCT ELISA Kit

Reactivity Human

Range 156 pg/ml - 10,000 pg/ml

Sensitivity

Tested Applications ELISA

Immunogen Expression system for standard: NS0; Immunogen sequence: V29-L361


Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only.

Sample Types cell culture supernatants, serum and plasma (heparin, EDTA).

Uniprot ID [Q16769](#)

Expiration Date Please enquire.



Human QPCT PicoKine
ELISA Kit