


Product Datasheet

Human PCOLCE ELISA Kit orb610885

Description	Human PCOLCE ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Sensitivity	
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q26-D449
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only.
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q15113
Expiration Date	Please enquire.



Human PCOLCE PicoKine
ELISA
Kit