


Product Datasheet

Mouse PCOLCE ELISA Kit orb610886

Description	Mouse PCOLCE ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Sensitivity	
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q25-A468
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only.
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q61398
Expiration Date	Please enquire.



Mouse PCOLCE PicoKine
ELISA
Kit