


Product Datasheet

Mouse Plasma Kallikrein/KLKB1 ELISA Kit orb610887

Description	Mouse Plasma Kallikrein/KLKB1 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Sensitivity	
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G20-A638
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only.
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	P26262
Expiration Date	Please enquire.



Mouse KLKB1 PicoKine
ELISA Kit