

Product Datasheet

GIPC1 antibody orb612191

Description	GIPC1 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Predicted Reactivity	Bovine, Human, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human GIPC1 (171-270/333 aa)
Target	GIPC1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	O14908
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.