

Product Datasheet

IRF3 antibody orb612212

Description	IRF3 antibody
Species/Host	Mouse
Predicted Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human IRF3
Target	IRF3
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG1
Clonality	Monoclonal
MW	47 kDa
Clone ID	2C4
Uniprot ID	Q14653
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IF=1:20-50
Expiration Date	12 months from date of receipt.