

Product Datasheet

CPT2 antibody orb612213

Description	CPT2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human CPT2
Target	CPT2
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	71 kDa
Clone ID	3B10
Uniprot ID	P23786
Dilution Range	WB=1:500-1000, IHC-P=1:50-200, ICC=1:50, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.