



Product Datasheet

Huntingtin antibody orb612215

Description	Huntingtin antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Predicted Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Hazard Information	Non-Toxic
Immunogen	Recombinant protein within human Huntingtin (1-150/3142 aa)
Target	Huntingtin
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	347 kDa
Clone ID	1F10
Uniprot ID	P42858
Dilution Range	IHC-P=1:50-200, WB=1:200-500, IHC-F=1:50-200, ICC=1:20-100
Expiration Date	12 months from date of receipt.



IHC-P analysis of mouse cerebellum using Huntingtin Monoclonal Antibody at 1:200 dilution.



IHC-P analysis of rat brain using Huntingtin Monoclonal Antibody at 1:200 dilution.