



Product Datasheet

GLUT1 antibody orb612216


Description	GLUT1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human GLUT1
Target	GLUT1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	54 kDa
Clone ID	2D5
Uniprot ID	P11166
Dilution Range	IHC-P=1:50-200, WB=1:300-500, IHC-F=1:50-200, ICC=1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin-embedded human liver tissue using anti-Glucose Transporter GLUT1 Antibody.



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue using anti-Glucose Transporter GLUT1 Antibody.



IHC-P analysis of rat liver using GLUT1 Monoclonal Antibody at 1:200 dilution.