



Product Datasheet

VEGFR1 antibody orb612220

Description	VEGFR1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human VEGFR1
Target	VEGFR1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	147 kDa
Clone ID	3A6
Uniprot ID	P17948
Dilution Range	WB=1:1000-2000, IHC-P=1:50-200, IHC-F=1:50-200, ICC=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P analysis of human brain using VEGFR1 Monoclonal Antibody at 1:200 dilution.



Western blot: Lane 1: Placenta (Mouse); Lane 2: Cerebrum (Mouse); Lane 3: Raw264.7 (Mouse) Cell; Lane 4: U2os (Human) Cell; Primary: Anti-VEGFR1 at 1:1000 dilution.