



Product Datasheet

PKC Delta antibody orb612223


Description	PKC Delta antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human PKC delta
Target	PKC delta
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	77 kDa
Clone ID	2C4
Uniprot ID	Q05655
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P analysis of mouse lymphoid using PKC delta Monoclonal Antibody at 1:200 dilution.



IHC-P analysis of rat cerebellum using PKC delta Monoclonal Antibody at 1:200 dilution.



Western blot: Lane 1: HeLa (Human) Cell; Lane 2: Jurkat (Human) Cell; Lane 3: MCF-7 (Human) Cell; Primary: Anti-PKC delta at 1:1000 dilution.