



Product Datasheet

S100 Beta antibody orb612224

Description	S100 Beta antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human S100B
Target	S100 beta
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	10 kDa
Clone ID	26C3
Uniprot ID	P04271
Dilution Range	WB=1:500-2000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P analysis of mouse brain using S100 beta Monoclonal Antibody at 1:50 dilution.



IHC-P analysis of rat brain using S100 beta Monoclonal Antibody at 1:200 dilution.



Western blot: Lane 1: zebrafish tissue; Primary: Anti-S100 beta at 1:500 dilution.