



Product Datasheet

Tubulin Beta-III antibody orb612226

Description	Tubulin Beta-III antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human Tubulin beta-III
Target	Tubulin beta-III
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	50 kDa
Clone ID	3B7
Uniprot ID	Q13509
Dilution Range	WB=1:1000-2000, IHC-P=1:50-200, ICC=1:50, Flow-Cyt=1:50, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: PC-12 Cell; Lane 2: HeLa Cell; Lane 3: SH-SY5Y Cell; Primary: Anti-Tubulin beta-III at 1:500 dilution.



Western blot: Lane 1: Mouse Cerebrum; Lane 2: Mouse Cerebellum; Lane 3: Rat Cerebrum; Lane 4: Rat Cerebellum; Lane 5: Human SH-SY5Y cell; Lane 6: Human U251 cell; Lane 7: Human U87MG cell; Primary: Anti-Tubulin beta-III at 1:1000 dilution.