



Product Datasheet

Bcl-2 antibody orb612229

Description	Bcl-2 antibody
Species/Host	Rabbit
Reactivity	Human
Predicted Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human Bcl-2
Target	Bcl-2
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	26 kDa
Clone ID	5F7
Uniprot ID	P49950
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/Test, IHC-P=1:100-500, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



Flow cytometry of THP-1. Primary Antibody (green line): Rabbit Anti-Bcl-2 antibody; Dilution: 1ug/Test.



Western blot: Lane 1: Jurkat (Human) Cell; Primary: Anti-Bcl-2 at 1:1000 dilution.