


Product Datasheet

RAB27A antibody orb612233

Description	RAB27A antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human RAB27A
Target	RAB27A
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG1
Clonality	Monoclonal
MW	25 kDa
Clone ID	9C6
Uniprot ID	P51159
Dilution Range	WB=1:500-2000, Flow-Cyt=1:20-50
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: K562 (Human) Cell; Lane 2: Jurkat (Human) Cell; Primary: Anti-RAB27A at 1:1000 dilution.