

Product Datasheet

Beta Amyloid 1-16 antibody orb612240

Description	Beta Amyloid 1-16 antibody
Species/Host	Rabbit
Predicted Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from rat beta-Amyloid (1-16/42 aa)
Target	beta Amyloid 1-16
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	4.3 kDa
Uniprot ID	P05067
Dilution Range	ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.