


Product Datasheet

ADAM28 antibody orb612241

Description	ADAM28 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human ADAM28 (521-620/775 aa)
Target	ADAM28
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	87 kDa
Uniprot ID	Q9UKQ2
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of human gastric using ADAM28 Polyclonal Antibody at 1:400 dilution.