


Product Datasheet

YAP1 antibody orb612243

Description	YAP1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Gallus, Porcine, Rabbit, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human YAP1 (51-150/504 aa)
Target	YAP1
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	55 kDa
Uniprot ID	P46937
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:293T(Human)
Cell;Testis
(Mouse);Primary:Anti-YAP1
at
1:1000 dilution.