

## Product Datasheet

### ZAR1 antibody orb612250

|                      |                                                                                  |
|----------------------|----------------------------------------------------------------------------------|
| Description          | ZAR1 antibody                                                                    |
| Species/Host         | Rabbit                                                                           |
| Reactivity           | Mouse                                                                            |
| Predicted Reactivity | Bovine, Human, Porcine, Rabbit, Rat, Sheep                                       |
| Conjugation          | Unconjugated                                                                     |
| Tested Applications  | WB                                                                               |
| Hazard Information   | Non-Toxic                                                                        |
| Immunogen            | KLH conjugated synthetic peptide derived from human ZAR1 (261-360/424 aa)        |
| Target               | ZAR1                                                                             |
| Form/Appearance      | Liquid                                                                           |
| Concentration        | 1mg/ml                                                                           |
| Preservatives        | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.            |
| Storage              | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note                 | For research use only.                                                           |
| Isotype              | IgG                                                                              |
| Clonality            | Polyclonal                                                                       |
| MW                   | 46 kDa                                                                           |
| Uniprot ID           | <a href="#">Q86SH2</a>                                                           |
| Dilution Range       | WB=1:500-2000                                                                    |
| Expiration Date      | 12 months from date of receipt.                                                  |