


Product Datasheet

ALKBH5 antibody orb612254

Description	ALKBH5 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Predicted Reactivity	Bovine, Canine, Human, Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human ALKBH5 (1-100/394 aa)
Target	ALKBH5
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	Q6P6C2
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot: Testis (Rat); Lung (Mouse); Primary: Anti-ALKBH5 at 1:1000 dilution.