

Product Datasheet

CD62P antibody orb612255

Description	CD62P antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human sp-selectin (501-600/830 aa)
Target	CD62P
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	86 kDa
Uniprot ID	Q01102
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.