



Product Datasheet

C22orf32 antibody orb612257

Description	C22orf32 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Predicted Reactivity	Bovine, Equine, Human, Porcine, Rabbit, Sheep
Conjugation	Unconjugated
Tested Applications	IHC-P
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human C22orf32 (41-107/107 aa)
Target	C22orf32
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	affinity purified by Protein A
MW	7kDa
Uniprot ID	Q9H4I9
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P image of rat pancreas tissue using C22orf32 antibody



IHC-P image of mouse pancreas tissue using C22orf32 antibody