



Product Datasheet

TBX3 antibody orb612259

Description	TBX3 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Predicted Reactivity	Human, Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-P
Hazard Information	Non-Toxic
Immunogen	KLH conjugated synthetic peptide derived from human TBX3 (181-280/74 aa)
Target	TBX3
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	82 kDa
Uniprot ID	O15119
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of mouse kidney using TBX3 Polyclonal Antibody at 1:200 dilution.



IHC-P analysis of rat liver using TBX3 Polyclonal Antibody at 1:200 dilution.