




Product Datasheet

TIMP1 antibody orb29164

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	TIMP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23171
Clone ID	RB14602
Uniprot ID	P01033
NCBI	7076, NP_003245.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, IF/ICC: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin-fixed paraffin embedded human pancreas tissue using TIMP1 antibody (C-term)



Confocal immunofluorescence analysis of A2058 cell using TIMP1 antibody (C-term)